

opto

Science in Sight



NEW
Flex
max

System for surgical
microscope for using with
wide field non-contact lens

MAXIMUM FLEXIBILITY IN RETINA AND VITREOUS SURGERY



NEW
Flex
max


MAXIMUM FLEX/8/LITY

The new Opto Flex Max comes with the new Opto SR (Super Range) wide field lens with high quality vision of the peripheral retina almost to the ora serrata. This new lens with the Opto Inverter is the ideal tool for vitreo retina surgeries.


The HM Lens is for vision of the macula. It provides a high magnification vision with a good deep of focus.

The new Opto Flex Max also provides comfort and freedom for the surgeon, and extreme fast procedures.


- Easy for handling, install and remove from microscope.
- Compatible with most microscopes.
- Surgeries with noassistance to hold lens.
- Sterilizable in autoclave.



Quick and easy to tilt



Eye protecting system



ilt turns 360°



Comes inside a carrying case with variable rings for different microscopes

* Opto manufactures Zeiss and Leica version
SPECIFICATIONS SUBJECTTO CHANGE WITHOUT PRIOR NOTICE OPTO FLEXMAX VITRECTOMY SYSTEM - ANVISA # 10309579006 Copyright• 2022 Opto Eletronica SA All rights reserved

OPTIONAL ITEMS



Opto Image Inverter* (comes with 810nm laser safety filter)



Opto green laser 532nm safety filter for microscopes*



Autoclavable Sterilization/Disinfection Case



MASTER GLOBAL DISTRIBUTOR

Contact: Juan Carlos Torres
info@jkmedicalinc.com
c: +1(407)342-8575
www.jkmedicalinc.com

